DENTSPLY

ECLIPSE®

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**Laboratory Directions For Use
Processing Eclipse® Prosthetics Resins in the
Austenal® Enterra™ VLC Curing Unit**

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EC REP
DeguDent GmbH
P.O. Box 1364 • 63403 Hanau
GERMANY
Tel. +49/6181/5950

© 2006 DENTSPLY International, Inc. All rights reserved.

Printed in U.S.A. 905577 Rev. C 8/06

CE

R_Xonly

Indications for Use: Eclipse® Prosthetic Resins are indicated for full baseplates, occlusal splints, nightguards, and provisional partial dentures when using the Enterra™ VLC Curing Unit.

Contraindications: See precautions and warnings.

Warnings: This product contains polymerizable materials. It does not contain methyl methacrylate (MMA) monomer. It is possible that some dental health care workers may develop an allergy or sensitivity to the product. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice. Do not immerse Eclipse appliances in boiling water.

Eclipse Bonding Agent should be used in accordance with instructions to prevent debonding.

Precautions:

1. This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
2. Users with special skin problems, cuts or abrasions may choose to wear protective gloves.
3. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
4. Since the Eclipse resins will become the final prosthesis, cleanliness is very important to an acceptable outcome. Wash hands before using uncured Eclipse materials.
5. Store at 60-80° F (16-27° C), away from moisture and in a light-safe container. The materials may polymerize prematurely if stored at excessively high temperature or if exposed to a visible or UV light source. Do not expose to direct sunlight.

6. Unused baseplate and clear resin materials should be cured prior to disposal.
7. Items removed from the curing unit are hot! Use heat-resistant gloves.
8. Do not remove prostheses from models until cool.
9. If the resin materials will not release easily from the package material (may occur in warm ambient conditions), place the package in a refrigerator/freezer for a few seconds to facilitate removal.
10. The following disinfectants have been used successfully with cured and uncured Eclipse resins:
 - Glutaraldehyde-based disinfectants
 - Quaternary ammonium chlorides
 - Quaternary ammonium chlorides/isopropyl alcohol (low concentration)
 - Spray-based phenolics

Use of other disinfectant products may not be suitable with Eclipse. Avoid disinfectants containing iodine.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding these materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

NOTE: Stop and read operating manuals for equipment before continuing to read Laboratory Instructions.

Instructions: Nightguard and Occlusal Splint:

1. Block out any undercuts on the model that could potentially lock the appliance in the mouth. A

duplicate model should be made. Normally the model will be damaged or destroyed during removal of the appliance after curing.

2. Place model in the Austenal® Enterra VLC Curing unit making sure the top of the model is below the maximum height line (row of holes in the back wall of the chamber). If any part of the model is above the maximum height line, remove stone from bottom of the model until the entire model is below the maximum height line.
3. Place a temperature indicator on the heel of the dry model. Paint a thin layer of Al-Cote® Separating Agent on the surface of the dry model and allow separator to dry completely.
4. Warm the model using heat lamps or in an oven set at 131°F (55°C). Heating time will vary depending on the number of models in the oven and the water content of the gypsum. Wait until the Temperature Indicator (provided) turns black on the heel of the model. This indicates when model is at correct temperature.
5. Place the rounded side of the Clear Baseplate Resin arch on the incisal/occlusal surface of the model and carefully adapt the resin to the desired shape. Trim any excess material with DENTSPLY Ney® Wax Pencil Pro (temperature setting at 3 o'clock position is suggested).

Note: Do not extend the nightguard below the height of contour as this could result in the final appliance locking into the mouth and difficulty in removal.
6. Before placing back on the articulator to check occlusion, coat the opposing dentition with Eclipse Model Release Agent. Then place the model on the articulator and gently tap into occlusion until the pin is down. Using the incisal guide pin and table, the nightguard can be put through side-to-

side excursions if desired. Remove model with nightguard for processing.

7. Apply Eclipse Air Barrier Coating to the entire surface of the nightguard, and place it in the Austenal Enterra VLC Curing Unit within 10 minutes; then start the nightguard cure sequence (Eclipse Family – Nightguard Cure Cycle). Be sure to center the device on the turntable.
8. When the curing cycle is complete, remove the nightguard and model from the Enterra VLC Curing Unit and allow to cool to room temperature.
9. Soak the nightguard and model in tap water to rehydrate model to facilitate nightguard removal from the model.
10. Finish and polish in the usual manner.

Note: The appliance should fit passively on the master model.

Heat & Seat™ (Compliant) Splint/Nightguard:

1. Place model in the Austenal Enterra VLC Curing unit to ensure the top of the model is below the maximum height line (row of holes in the back wall of the chamber). If any part of the model is above the maximum height line, remove stone from bottom of the model until the entire model is below the maximum height line.
2. Dry the model completely and block out any severe undercuts on lingual, buccal and interproximal areas and make a duplicate model prior to fabrication when necessary for processing.
3. Paint the model with Al-Cote® Separating Agent and allow it to dry. Separator **MUST** be dry before fabricating.
4. Place the model in an oven to warm slightly OR warm model with a hot air gun for about 1 minute. Model must be warm for

5. Remove Eclipse Resilient material from package by first holding the package under cold running water for 15-30 seconds to facilitate easy removal. Material can be refrigerated.
6. On the warm model, adapt the Eclipse Resilient Resin to the teeth and occlusal surface based on the design of the appliance. You can allow cusps to penetrate the Resilient Resin depending on design. If the model has cooled down, warm it with a hot air gun and continue to adapt the Eclipse Resilient Resin. Eclipse Resilient Resin usually will not need to be placed below the height of contour of the teeth as with other soft materials. Any excess material can be removed using a sharp instrument.
7. Once the Resilient material is in position on the occlusal surface, remove Clear Eclipse Baseplate Resin from package and cut in half lengthwise. You normally do not need a complete arch of material for the occlusal portion of the compliant nightguard. Adapt the Clear Eclipse Baseplate Resin directly on top of the resilient material already in place making sure not to displace the Resilient resin. If material and model have cooled down, warm in an oven or with a hot air gun to facilitate easier adaptation.
8. Before placing back on the articulator to check occlusion, coat the opposing dentition with Eclipse Model Release Agent. Then place the model on the articulator and gently tap into occlusion until the pin is down. Using the incisal guide pin and table, the nightguard can be put through side-to-side excursions if desired.
9. Smooth the junction of the Clear Eclipse Baseplate and the Resilient Resins using the DENTSPLY Ney Wax Pencil Pro (temperature setting at 3 o'clock position is suggested). Further smoothing can be accomplished with a hot air gun as well.
10. Remove the model from the articulator and paint the appliance with Eclipse Air Barrier Coating. Immediately process in the Austenal Enterra VLC Curing Unit, using Eclipse Family-Nightguard Cure Cycle. Be sure to center the device on the turntable.
11. Upon completion of the curing cycle, allow to cool to room temperature before rinsing off the Eclipse Air Barrier Coating. Remount on the articulator to verify/adjust occlusion as necessary.
12. Soak splint and model in warm water for 5 minutes to facilitate removal from the model. This warming allows the Eclipse Resilient material to become slightly flexible for easy removal from the model.
13. Wash off Eclipse Air Barrier Coating with water and a soft brush.
14. Finish and polish. **Note:** For seating and removal of Eclipse Resilient appliances it is necessary to warm the appliance with warm water. This warming allows the Eclipse Resilient material to become slightly flexible for easy removal or insertion.

Optional Two-cure Method

Follow steps 1-6, but after the resilient material is applied process the Eclipse Resilient Resin in the Austenal Enterra VLC Curing Unit, using Eclipse Family-Resilient cure cycle. **NOTE:** Be sure to center the device on the turntable. Do NOT apply Eclipse Air Barrier Coating to the resilient material before processing. After curing remove from Austenal Enterra VLC Curing Unit and continue with steps 7-14 above.

Stabilized Baseplate for Eclipse/Lucitone 199® Hybrid Denture:

1. Box and pour the master model. The land area of the model should be 1mm or less above the depth of the sulcus to allow light penetration for complete curing.
2. Carve the post-dam in the master model if not included in the impression.
3. Place model in the Austenal Enterra VLC Curing unit to ensure the top of the model is below the maximum height line. If any part of the model is above the maximum height line, remove stone from bottom of the model until the entire model is below the maximum height line.
4. Due to undercuts, a duplicate or working model should be made. The baseplate must be cured on the master model. Normally the master model will be damaged or destroyed during removal of the baseplate after curing.
5. Place a temperature indicator on the heel of the dry model. Paint a thin layer of Al-Cote® Separating Agent on the surface of the dry master model and allow separator to dry completely. Note: Since the baseplate may become part of the final prosthesis, cleanliness is very important for an acceptable outcome.
6. Warm the master model to 120-130°F (49-54°C). This can be accomplished under heat lamps or in an oven set at 131°F (55°C). Heating time will vary depending on the number of models in the oven and the water content of the gypsum. Wait until the Temperature Indicator (provided) turns black on the heel of the model. This indicates when model is at correct temperature.
7. Place the rounded side of the upper or lower arch of Eclipse Baseplate Resin onto the ridge and allow the Baseplate Resin to warm on the model for 30 seconds. Carefully adapt the labial/buccal surfaces to avoid air entrapment. Begin adaptation at the crest of the ridge and slowly and carefully, using your thumbs and fingers, adapt the material down the buccal slope and into the vestibule last. Follow the same procedure for the palatal/lingual section. For the maxillary baseplate, carefully join the two sides at the midline without trapping air.
8. Once the baseplate is adapted to the model, apply Eclipse Air Barrier Coating to the entire surface of the baseplate, and place it in the Austenal Enterra VLC Curing Unit within 10 minutes; then start the baseplate cure sequence (Eclipse Family – Flipper/Base Cure Cycle). Be sure to center the device on the turntable.
9. When the curing cycle is complete, remove the baseplate/model from the Austenal Enterra VLC Curing Unit and allow the baseplate/model to cool to room temperature before proceeding.
10. Soak the baseplate and master model in tap water to rehydrate model to facilitate baseplate removal from the master model. If undercuts are present, section the master model carefully with a separating disc from the underside prior to removal.
11. If tissue side shade is unacceptable (slightly orange in color), apply Eclipse Air Barrier Coating to the tissue side of the baseplate. Place baseplate in the Austenal Enterra VLC Curing Unit tissue side up. Start the cure sequence (Eclipse Family – BP Tissue Side Cure Cycle). Be sure to center the device on the turntable.
12. When the curing cycle is complete, remove the baseplate from the Austenal Enterra VLC curing unit and allow the baseplate to cool on the bench to room temperature.
13. Wash Eclipse Air Barrier Coating off of baseplate using water and a brush.

USA GB

14. Trim and smooth the baseplate borders. If the master model is damaged during removal of the baseplate, the duplicate model may require modification to allow seating of this baseplate for articulation.
- Refer to the Eclipse Wax Remover Directions for Use to fabricate a hybrid full denture using the processed Eclipse baseplate and Lucitone 199 denture base resin.

PROVISIONAL PARTIAL DENTURES:

1. Eclipse Bonding Agent needs to be used to prepare the teeth for fabrication in a provisional appliance. Follow the bonding agent instructions for use for tooth preparation.
2. When using porcelain teeth with pins and/or vent holes, no tooth preparation is required but you will need to flow baseplate resin into the vent holes and around the pins prior to tooth placement.
3. Be sure that entire model is below the maximum height line.
4. If wire clasps are utilized, position them onto the warm model with a small amount of Eclipse baseplate resin, underneath each retention loop or tail and glue the buccal clasp arms in place on the model.
5. With teeth prepared using Eclipse Bonding Agent and clasps in position if required, apply a sufficient amount of the Eclipse Baseplate Resin to the warm model.
6. Adapt the resin to the warm model being careful not to create or trap any air bubbles. Smooth the material and set the teeth into the resin and contour.
7. Once appliance is shaped and smoothed paint it with Eclipse Air Barrier Coating and process using Eclipse Family – Flipper/Base Cure Cycle. Be sure to center the device on the turntable.
8. When the curing cycle is complete, remove the appliance from the Austenal Enterra VLC Curing Unit and allow to cool to room temperature.
9. Soak the appliance and master model in tap water to rehydrate model to facilitate removal from the master model.
10. If tissue side is unacceptable (slightly orange in color), apply Eclipse Air Barrier Coating to the tissue side of the appliance. Place appliance in the Austenal Enterra VLC Curing Unit tissue side up. Start the cure sequence (Eclipse Family BP Tissue Side Cure Cycle). Be sure to center the device on the turntable.
11. When the curing cycle is complete, remove the appliance from the Austenal Enterra VLC Curing Unit and allow the appliance to cool on the bench to room temperature.
12. Wash Eclipse Air Barrier Coating off of the appliance using water and a brush.
13. Finish and polish in the usual manner.

5

USA GB

- 3 o'clock position is suggested) to flow the Eclipse Baseplate Resin onto the prepared area.
- Carefully flow from one side to the other to avoid air entrapment.
 - It is vital that a melt interface occurs between the Baseplate Resin and the prepared surface.
 - Fill to contour or slightly below. Do not overfill; otherwise fit may be affected!
 - If multiple voids need to be repaired, locally warm the repaired area using a hot air gun to ensure that a melt interface occurs.
5. Apply Eclipse Air Barrier Coating to any exposed surfaces.
 6. Cure in the Austenal Enterra VLC Curing unit using the Eclipse
 7. Family - BP Repair cure cycle. Be sure to center the device on the turntable.
 7. Immediately remove the repaired baseplate from the Austenal Enterra VLC Curing Unit and place the repaired baseplate back on the model and allow to cool to room temperature.
 8. Remove Eclipse Air Barrier Coating with water and a brush.
 9. Finish the repair and polish in the usual manner.

Repair of Minor Voids in Baseplate:

1. For small voids, prepare the area with a coarse bur or diamond to enhance bonding. Do not leave feathered edges.
2. Thoroughly clean prepared surface as described below:
 - Wash prepared area with a brush and tap water and blow dry with clean compressed air.
 - Use brush and isopropyl or ethyl/denatured alcohol to clean the prepared area and allow to dry.
 - Make sure no dust remains on the surface.
3. Locally warm the prepared area with a hot air gun.
4. Use the DENTSPLY Ney Wax Pencil Pro (temperature setting at

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SECTION IV - FIRE AND EXPLOSION HAZARD DATA

- 4.1 Flash Point (Method Used) > 220° C (428° F)
- 4.2 Flammable (Explosive) Limits
LEL: NA UEL: NA
- 4.3 Extinguishing Media Dry chemical, CO₂, foam, water spray, and water fog
- 4.4 Special Fire Fighting Procedures and Equipment: Fireman should wear protective clothing and self-contained respirators if large amounts of material are involved in a fire.
- 4.5 Unusual Fire and Explosion Hazards: Irritating gases may be emitted when combusted.

Eclipse Denture Base Resin [Baseplate, Setup and Contour Resin], cont'd

SECTION V - HEALTH HAZARD DATA

- 5.1 Effects of Acute and Chronic over Exposure: Skin irritation is possible in individuals sensitive to (meth)acrylates and contact dermatitis (skin sensitization) can develop in some individuals.
- 5.2 Emergency and First Aid Procedures: Eyes: If material contacts eyes, flush eyes immediately with copious amounts of tap water for a minimum of 15 minutes and get medical attention. Skin: If skin irritation occurs, wash irritated area with soap and water. Remove clothing if it becomes contaminated with resin. Sensitization: To minimize chance of becoming sensitized, wash hands often when using material. Particularly sensitive individuals should avoid skin contact.

SECTION VI - REACTIVITY DATA

- 6.1 Unstable: Stable: X
- 6.2 Conditions to avoid (Stability): Material can polymerize if exposed to high temperatures or to direct sunlight.
- 6.3 Incompatibility (Materials to avoid): Oxidizing agents.
- 6.4 Hazardous Decomposition Products: Acid fumes may result if large amounts are combusted during a fire.
- 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X
- 6.6 Conditions to avoid polymerization: Direct sunlight, fluorescent lighting, high temperatures for extended time period.

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 Handling and Storage Precautions: This material should be handled and used only in a manner consistent with the product instructions. Do not store at high temperature or in sunlight. Store away from food.
- 7.2 Other Precautions: Visible light can cause polymerization.
- 7.3 Steps to be taken in case material is released or spilled: No Special procedure required.
- 7.4 Waste Disposal Procedures: Polymerize baseplate and clear resin by exposing material to sunlight until solid, then dispose of as a non-hazardous waste. Small (gram / oz.) quantities of set up resin and contour resin can be disposed of as non-hazardous waste.
- 7.5 Neutralizing Agent: NA

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 Respiratory Protection: Use a dust mask while grinding the cured material to avoid breathing dust.
- 8.2 Ventilation: Use local exhaust when grinding cured material.
- 8.3 Protective Gloves: For users with skin problems, cuts or abrasions - wear disposable protective gloves.
- 8.4 Eye Protection: Use safety glasses while grinding.
- 8.5 Other Personal Protective Equipment: None required

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